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Software for recording and evaluating the condition of members of the army and first responders

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The result is a specialized acquisition and analysis application designed for the collection, integration, and evaluation of multimodal time-series data related to movement activity and physiological responses of members of the army and first responders. The key principle is a unified data layer based on InfluxDB, which consolidates data streams from various sources (image data, inertial sensors, ECG, stabilometry) into a common timeline. The application enables storage, retrieval, export, and reproducible evaluation of data. It supports a real-time visualization mode for operational integrity checks during

measurement and an offline analysis mode for post hoc inspection and validation. In the camera-based branch, the MediaPipe Pose model is integrated for estimating body posture. Each measurement results in a database-stored session identified by a session_id tag, serving as a unified source for CSV export as well as for internal offline analysis within the application. The software provides a consolidated platform for the collection, monitoring, and evaluation of heterogeneous data streams, enabling the creation of validated datasets and supporting decision-making regarding users' physical fitness and health status.

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multimodal measurement, time series, stabilometry, posture, data analysis

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Abstract (EN)

The result is a specialized acquisition and analysis application for multimodal data collection and offline analysis of time-series data related to movement activity. The key principle is a unified data layer based on InfluxDB, which consolidates asynchronous data streams into a common timeline and enables auditable storage, repeated querying, export, and reproducible evaluation. The application is implemented in Python 3.10 with a Qt (PySide6) graphical user interface and a multithreaded architecture. In the vision branch, the MediaPipe Pose model is integrated for extracting body keypoints, and optional 3D triangulation is supported following stereo calibration. In real-time mode, visualization primarily serves for operational integrity checks of the measurement. The final output is a database-stored session identified by session_id, which can be used for CSV export as well as for internal offline analysis within the application.

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Keywords (EN)

multimodal measurement; time series; InfluxDB; real-time visualization; offline analysis; MediaPipe Pose; 3D triangulation; CoP; ECG; IMU

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
A screenshot of a Windows command prompt window. The title bar shows the path 'C:\Windows\System32\cmd.exe'. The window content displays the following text:

```
Microsoft Windows [Version 10.0.26100.7462]
(c) Microsoft Corporation. Všechna práva vyhrazena.

C:\Program Files\Influxd>influxd|
```

```
C:\Windows\System32\cmd.exe
\volfp\...\influxdbv2\engine\data\7c95a42434a2a68a\autogen\28\000000004-000000001.tsm"}
2026-01-15T10:37:07.704679Z info Beginning compaction {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "tsml_files_n": 3}
2026-01-15T10:37:07.705277Z info Compacting file {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "tsml_index": 0, "tsml_file": "C:\\Users\\..."}
2026-01-15T10:37:07.705277Z info Compacting file {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "tsml_index": 1, "tsml_file": "C:\\Users\\..."}
2026-01-15T10:37:07.705277Z info Compacting file {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "tsml_index": 2, "tsml_file": "C:\\Users\\..."}
2026-01-15T10:37:07.762127Z info Compacted file {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "tsml_index": 0, "tsml_file": "C:\\Users\\..."}
2026-01-15T10:37:07.762127Z info Finished compacting files {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "tsml_files_n": 1}
2026-01-15T10:37:07.762127Z info TSM compaction (end) {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "op_event": "end", "op_elapsed": "57.447ms"}
2026-01-15T10:37:07.879824Z info Compacted file {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "tsml_index": 0, "tsml_file": "C:\\Users\\..."}
2026-01-15T10:37:07.879824Z info Finished compacting files {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "tsml_files_n": 1}
2026-01-15T10:37:07.879824Z info TSM compaction (end) {"log_id": "10SJXPKG000", "service": "storage-engine", "engine": "tsml", "tsml_strategy": "full", "tsml_optimize": false, "op_name": "tsml_compact_group", "op_event": "end", "op_elapsed": "175.144ms"}
```

The screenshot displays the FluxusDB web interface. The 'DEVICE CONTROL' section at the top allows selecting between 'Movesense IMU+ECG' and 'FLIR Stereo System'. The 'PLATFORM' section shows the system is 'Ready' and provides fields for COM port (COM3), width (1000.00mm), and length (1000.00mm). Below this, the 'DB' section has a 'CONFIG' button highlighted with a red arrow. The 'RECORDING' section shows a subject ID of 'subject_01' and a recording status of 'REC' with a duration of '00:00.0'. On the right, a 'FluxusDB' configuration dialog box is open, showing fields for URL ('http://localhost:8086/'), Token (masked with dots), Org ('cvut'), Bucket ('MovArm'), and 'Save'/'Cancel' buttons.

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